

InvA colony PCR

Primers:

The primer used here came from:

Kim, J. S., G. G. Lee, J. S. Park, Y. H. Jung, H. S. Kwak, S. B. Kim, Y. S. Nam, and S.-T. Kwon. 2007. A novel multiplex PCR assay for rapid and simultaneous detection of five pathogenic bacteria: *Escherichia coli* O157:H7, *Salmonella*, *Staphylococcus aureus*, *Listeria monocytogenes*, and *Vibrio parahaemolyticus*. J. Food Prot. 70:1656–1662.

The sequences are (5' => 3'):

invA-F: GAATCCTCAGTTTTTCAACGTTTC

invA-R: TAGCCGTAACAACCAATACAAATG

Master Mix:

PCR Reagents [Concentration]	Vol. (μL) for 1x 50μL rxn	
dH ₂ O	31	
5x Go Taq Flexi buffer	10	
MgCl ₂ [25mM]	3	
dNTPs [10mM]	1	
primer: invA-F [12.5μM]	2	
primer: invA-R [12.5μM]	2	
Go Taq DNA polymerase	0.25	
TOTAL:	49.25	

Remember:

-Use GoTaq instead of Amplitaq

-refer to Colony PCR Protocol (#8.1.1.1.4) for DNA template preparation -add 1μL (or 1.6μL if using the multichannel pipettor) of dirty lysate to each corresponding master mix tube

-amplicon should be 678 bp

PCR Conditions:

94°C for 2 min [1X]

94°C for 30 s

60°C for 30 s 72°C [20X]

for 30 s

72°C for 5 min [1X]

Troubleshooting:

If colony PCR doesn't work, try making dirty lysates:

- Pick a colony with a toothpick and scrape into PCR tube, add 100µL of dH₂O
- Microwave dirty lysates for 30s to lyse cells